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

 ITFIT	SAMSUNG Lifestyle																																																
	Before Use <p>Thank you for using Rechargeable Electric Toothbrush. Before operating this product, please read the User's Manual completely and save them for future use.</p>																																																
	Feature of Product <ul style="list-style-type: none">• 5 Brushing Modes, store the mode lately used• 48,000 VPM• Use up to 180min. when fully charged• 6 Brush Heads• USB Rechargeable• IPX7 Waterproof																																																
	Product Information <p>Product Name: ITFIT Rechargeable Electric Toothbrush Model No. ITFITEX37 Material: ABS, Cylindrical Lithium-ion Battery, Electronic Parts Input (Charging Station): 5V/1A Power Consumption (Charging Station): 1.8W Battery Capacity: 800mAh (2.96Wh) / 3.7V Product Dimensions (Brush Head and Handle): 233x28x28mm Product Weight (Brush Head and Handle): Approx. 101g Travel Case Dimensions: 201x71x37mm Travel Case Weight: 83g Operating/Storage Temperature: 0°C-45°C Electric Toothbrush Waterproof Rating: IPX7 Country of Manufacture: China Manufacturing Date: 2023.12</p>																																																
Rechargeable Electric Toothbrush User Manual																																																	
	Product Appearance Introduction 																																																
	Mode Descriptions <table border="1"><thead><tr><th>Mode Names</th><th>Indicator Light</th><th>VPM (Vibrations Per Minute)</th><th>Vibration Frequency</th></tr></thead><tbody><tr><td>Clean</td><td>Clean</td><td>Approx. 35000</td><td>295±20Hz</td></tr><tr><td>White</td><td>White</td><td>Approx. 41000</td><td>342±24Hz</td></tr><tr><td>Polish</td><td>Polish</td><td>Approx. 48000</td><td>400±28Hz</td></tr><tr><td>Massage</td><td>Massage</td><td>Approx. 41000</td><td>342±24Hz</td></tr><tr><td>Sensitive</td><td>Sensitive</td><td>Approx. 35000</td><td>295±20Hz</td></tr></tbody></table> <p>5. Briefly press the power button for manual shutdown (automatic shutdown after 2 minutes). Notes:<ul style="list-style-type: none">• Within the initial 5 seconds before the mode starts, a short press of the power button shifts to the next mode. Beyond 5 seconds into the mode, a short press of the power button powers off the device.• Clean/Massage/Sensitive modes share the same VPM and frequency but vary slightly in vibration intensity and pattern.• The toothbrush features a built-in timer with a prompt every 30 seconds, reminding you to clean each quadrant of your mouth for approximately a quarter of the recommended brushing time.</p>	Mode Names	Indicator Light	VPM (Vibrations Per Minute)	Vibration Frequency	Clean	Clean	Approx. 35000	295±20Hz	White	White	Approx. 41000	342±24Hz	Polish	Polish	Approx. 48000	400±28Hz	Massage	Massage	Approx. 41000	342±24Hz	Sensitive	Sensitive	Approx. 35000	295±20Hz																								
Mode Names	Indicator Light	VPM (Vibrations Per Minute)	Vibration Frequency																																														
Clean	Clean	Approx. 35000	295±20Hz																																														
White	White	Approx. 41000	342±24Hz																																														
Polish	Polish	Approx. 48000	400±28Hz																																														
Massage	Massage	Approx. 41000	342±24Hz																																														
Sensitive	Sensitive	Approx. 35000	295±20Hz																																														
	Electric Toothbrush Usage Tips <ul style="list-style-type: none">• Start brushing from the edges of the teeth, gently guiding the electric toothbrush across the entire tooth.• Avoid applying excessive pressure with the brush head to reduce the risk of tooth sensitivity.• It's recommended to brush the lower teeth first, then the upper teeth. Brush the outer surfaces of the teeth before moving to the inner surfaces.• Brush 2-3 teeth at a time during the brushing process.• For a better use of the product, you may consult with a dentist.																																																
	Charging Indicator Light Guide <p>Placing the toothbrush on the charging station initiates the charging process. When charging, the charge indicator light flashes slowly;</p> <p>once fully charged, the charge indicator light remains solid.</p>																																																
	Warnings <ul style="list-style-type: none">• Consult a dentist before use if undergoing recent oral treatment.• Ensure the brush head is securely attached to the handle before each use; do not use without the brush head.• Do not use the product if unfamiliar with its operation to prevent potential harm.• Keep out of reach of children; not intended for use as a toy.• Cease use immediately if the product, charger, or charging device is damaged.• Avoid submerging the charger in water and refrain from prolonged exposure of the electric toothbrush to water to prevent malfunction, fire, electric shock, or explosion.																																																
	Notes: In order to extend battery life and ensure safer product usage, it is recommended to disconnect the power source after fully charging the product.																																																
	Daily Maintenance <ul style="list-style-type: none">• Keep the brush head clean to prevent the growth of bacteria and microbes.• After each use, rinse the brush head with water and pat it dry.• It is recommended to detach the brush head from the handle and rinse it at least once a week.• Regularly replace the brush head (recommended every 3 months). If you notice deformities or wear on the bristles, replace them promptly.• Rinse the brush head and handle with water. Use a soft cloth to wipe the charging station and travel case.																																																
	Eco-friendly Tips <p>This product contains a battery and is considered recyclable electronic waste. To protect the environment, please do not dispose of it with household garbage; instead, take it to the electronic waste recycling facility provided in your area.</p>																																																
	Declaration of Restricted Substances Content <table border="1"><thead><tr><th rowspan="2">Unit</th><th colspan="6">Restricted substances and its chemical symbols</th></tr><tr><th>Lead(Pb)</th><th>Mercury(Hg)</th><th>Cadmium(Cd)</th><th>Hexavalent Chromium(Cr6+)</th><th>Polybrominated Diphenyl Ethers (PBDEs)</th><th>Phthalates (Phthalic acid derivatives)</th></tr></thead><tbody><tr><td>PCB板</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr><tr><td>塑膠零件</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr><tr><td>電子零件</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr><tr><td>充電座</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr><tr><td>外觀材質</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td><td>○</td></tr></tbody></table> <p>Note 1: "Exceeding 0.1 wt%" and "exceeding 0.01 wt%" indicate that the percentage content of the restricted substance exceeds the reference percentage value of presence condition. Note 2: "○" indicates that the percentage content of the restricted substance does not exceed the percentage of reference value of presence condition. Note 3: "○" indicates that the restricted substance corresponds to the exemption.</p> 	Unit	Restricted substances and its chemical symbols						Lead(Pb)	Mercury(Hg)	Cadmium(Cd)	Hexavalent Chromium(Cr6+)	Polybrominated Diphenyl Ethers (PBDEs)	Phthalates (Phthalic acid derivatives)	PCB板	○	○	○	○	○	○	塑膠零件	○	○	○	○	○	○	電子零件	○	○	○	○	○	○	充電座	○	○	○	○	○	○	外觀材質	○	○	○	○	○	○
Unit	Restricted substances and its chemical symbols																																																
	Lead(Pb)	Mercury(Hg)	Cadmium(Cd)	Hexavalent Chromium(Cr6+)	Polybrominated Diphenyl Ethers (PBDEs)	Phthalates (Phthalic acid derivatives)																																											
PCB板	○	○	○	○	○	○																																											
塑膠零件	○	○	○	○	○	○																																											
電子零件	○	○	○	○	○	○																																											
充電座	○	○	○	○	○	○																																											
外觀材質	○	○	○	○	○	○																																											
	SAMSUNG C&T TAIWAN CORPORATION Address: 110 Rm. 2203 22F, No.333, Keelung Rd., Sec.1 Taipei, Taiwan TEL: +886 2 2728 8500 E-mail: itfittw@samsung.com																																																

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